Corning® Stripettor™ Ultra Pipet Controller

Instruction Manual

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Manual de Instrucciones

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





Table of contents

1. Introduction	2
2. Work safety instructions	2
3. Limitations of use	3
4. Switching on	3
5. Aspirating and dispensing liquids	3
6. Troubleshooting	5
7. Replacing the filter	6
8. Charging the batteries	6
9. Maintenance	7
10. Components	7
11. Ordering information	8
12. Spare parts	8
13. Product disposal	8

1. Introduction

The pipet controller is a device intended for pipetting liquids with the use of measuring pipets. It can work with all types of glass or plastic pipets in the volume range from 0.5 mL to 100 mL.

Two dispense modes permit selection of dispensing intensity depending on the user's needs. The selected setting of the pipet controller mode is shown on the display (Figure 1E).

Figure 1 shows the external parts of the pipet controller with a description of the materials used.

Display indicators are shown in Figure 2.

2. Work safety instructions

 **WARNING! Risk of injury**

CAUTION: Risk of damage to the device or errors in pipetting of liquids

Before starting the work with the pipet controller every user should read these operating instructions carefully.

CAUTION:

- Using the device inconsistently with the operating instructions may result in damaging the device.
- The device should be serviced only at an authorized service center, otherwise the manufacturer will be relieved from any liability under the warranty.
- Only original spare parts and accessories, recommended by the manufacturer, shall be used.
- Only the original charger, supplied by the manufacturer, shall be used for charging the batteries.
- In case of incorrect functioning of the pipet controller, work shall be stopped. The device shall be cleaned according to section 9 and sent for repair to an authorized service center.
- In the case of mechanical damage to the casing, the device shall be immediately sent for repair to an authorized service center.
- The use of excessive force during work shall be avoided.

 **WARNING!**

- During the work with the pipet controller general safety regulations regarding risks related with laboratory work should be observed. Protective clothing, goggles, and gloves should be worn.
- The pipet controller shall be used only for measuring liquids in conditions specified by the manufacturer, which are limited due to the chemical and mechanical resistance of the device, as well as the user safety.
- The information and instructions provided by the manufacturers of the reagents must be observed.

Note: The pipet controller is equipped with liquid vapors exhaustion system which protects against corrosion to ensure long instrument life.

3. Limitations of use

- The pipet controller shall not be used for measuring substances with vapours of which damage the following plastics: PP, SI, EPDM, POM.
- The pipet controller shall not be used in an environment where explosion risk is present.
- Flammable liquids shall not be measured – in particular substances with flash-point below 0°C (ether, acetone).
- The pipet controller shall not be used for drawing acids with a concentration above 1 mol/L.
- The pipet controller shall not be used for drawing solutions with a temperature above 50°C.
- The pipet controller may work in temperature range from +10°C to +35°C.

The pipet controller is suitable for general laboratory use only. It must be used only by personnel who know the health risks associated with the substances that are normally used with this instrument.

4. Switching on

The pipet controller is switched on by pressing any of the buttons (Figure 1A, B, C, D). The display will show the selected aspiration mode, dispense mode and the batteries level indicator. Examples of the display indications are shown in (Figure 5). The batteries are discharged and require charging if the indicator does not display "bars" (Figure 5.4). When the batteries are fully charged, the indicator displays three "bars" (Figure 5.1).

- The pipet controller switches off automatically if not used for 5 minutes.
- The pipet controller may be charged only with the original charger.
- The mains voltage shall conform with the specification on the charger.
- Charging shall be done in accordance with section 8 of the instruction manual.

5. Aspirating and dispensing liquids

Attaching a pipet

CAUTION: Before attaching a pipet, check whether the pipet is not damaged, has no dents or sharp edges in the gripping part. Check whether the gripping part is dry.

The pipet shall be gripped as close to the upper end as possible and carefully inserted into the pipet holder until resistance is noticed (Figure 3A).

WARNING!

Do not apply excessive force so as not to damage thin pipets and avoid risk of injury. A pipet that has been correctly attached and sealed in the holder should not tilt to the sides.

After attaching a pipet, hold the pipet controller in a vertical position. It is not recommended to leave the device with a pipet attached for a longer period, for example overnight or over a weekend.

CAUTION: Do not put aside the pipet controller if there is liquid in the pipet.

Filling the pipet

Before aspirating is started, set the speed by pressing the SPEED switch (Figure 1C) until the display shows the correct speed (Figure 1E).

- HIGH speed - fast aspirating (Figure 2.1),
- LOW speed - slow aspirating (Figure 2.2).

It is recommended to set the LOW speed when working with pipets of volumes up to 5 mL, and the HIGH speed for pipets of volumes greater than 5 mL. Holding the pipet controller in a vertical position, immerse the pipet end in the liquid to be drawn up (Figure 3B), and press the aspiration button gently. The speed depends on how deep the aspirating button has been pressed. The deeper the button is pressed the faster the liquid is aspirated into the pipet.

It is recommended to draw a slightly greater liquid volume than required (due to meniscus above the required volume mark), adjusting the aspiration speed, so as not to overfill the pipet.

Setting the volume

After the pipet is filled, dry the outside surface with absorbent paper that does not leave impurities. Then set the required liquid volume precisely. Pressing the dispense button gently (Figure 3C), dispense the excessive liquid from the pipet until the meniscus of the liquid aligns exactly with the required volume mark on the pipet.

Emptying the pipet

Holding the vessel in an inclined position, place the pipet end in contact with the vessel wall and press the dispense button gently (Figure 3C). The dispensing intensity may be adjusted depending on how deep the dispense button has been pressed. The deeper the button is pressed the faster the outflow of liquid from the pipet.

The pipet controller has two dispense modes. The dispense mode is selected by successive pressing the MODE switch (Figure 1D) until the display shows the correct mode (Figure 1E).

- Gravity mode marked with the  icon on the display (Figure 2.4) – dispensing is effected in gravity mode, which means that the liquid flows out of the pipet by its own weight.
- Blow out mode marked with the  icon on the display (Figure 2.5) – dispensing is effected in gravity mode, however, when the dispense button is pressed to the middle position, the pump is started and fast emptying of the pipet with a blow out is effected.

CAUTION:

During gravimetric dispensing the pipet is not completely emptied due to the characteristics of pipets used with the pipet controller.

6. Troubleshooting

If during your work the functioning of the pipet controller is incorrect, check the cause and remove the fault.

Problem	Possible Cause	Action
The pipet falls out (the holding force of the pipet is too small), or tilts to the side too much.	The pipet holder is dirty or wet (Figure 1G).	Take out the pipet holder, then clean, wash, and dry it.
	The pipet holder is damaged.	Replace the holder with a new one.
The pump is working, but the pipet controller aid does not draw liquid or draws liquid very slowly.	The filter is dirty (Figure 1H).	Take out the pipet holder, take out the filter, if it is dirty, replace it with a new one.
	The pipet holder and/or the connector gasket is damaged (Figure 1J).	Replace the mechanically damaged elements with new ones.
Liquid leaks from the pipet (the aspiration and the dispense buttons are not pressed).	The pipet is damaged.	Check the pipet for damage (cracks, dents). If present, replace the pipet with a new one.
	The pipet is inserted incorrectly.	Check whether the pipet has been correctly inserted in the holder.
	The pipet holder, the filter, or the connector gasket is installed incorrectly.	Check whether all parts are present and correctly installed.
	The pipet holder and/or the connector gasket is damaged (Figure 1G, Figure 1J).	Replace the mechanically damaged elements with new ones.

If the above actions do not help, the device shall be sent to the nearest authorized service center.

Before being serviced, the pipet controller should be cleaned and decontaminated. Written details including the precise specification of solutions used and the type of laboratory in which the device was used, should be sent with the product.

7. Replacing the filter

CAUTION:

The work safety instructions given in section 2 shall be observed when disassembling the pipet controller.

The filter replacement is necessary, if drawing efficiency deterioration is observed. The direct reason may be a dirty filter after a long period of use.

In order to replace the filter:

- Remove the pipet.
- Unscrew the nose piece (Figure 4A).
- Remove the membrane filter (Figure 4A) and the pipet holder (Figure 4B).
- Rinse the holder using a wash bottle (Figure 4C).
- Blow liquid out of the holder and set it aside until it is completely dry.
- Install new membrane filter (Figure 4D) and assemble the device in reverse order.

8. Charging the batteries

CAUTION!

The pipet controller may be charged only with the original charger. The main voltage must conform with the specification on the charger.

Using chargers other than the original one may damage the battery.

The pipet controller is delivered with 3 NiMH type AAA batteries.

The batteries can be replaced easily if necessary after the battery cover is taken off (Figure 1K). The manner in which the batteries are arranged is shown in Figure 8.

Charging:

1. Charging temperature: 10°C to 35°C.
2. Charging the battery is carried out through a charger (power supply) by direct connection to the main power. Batteries charging is indicated by successive lighting of "bars".
3. Full charging time: 7 to 8 hours.
4. The batteries are charged when all 3 "bars" are displayed simultaneously (Figure 5.1).

When the batteries are charged, the charging circuit disconnects automatically.

The service life of the batteries: approx. 1,000 charging cycles, if used correctly. It is not possible to overcharge the batteries if all instructions of the manufacturer are followed.

WARNING!

If disposable batteries are used, it is unacceptable to connect the device to the charger. In order to prolong the life span of the rechargeable batteries, the following rules should be followed:

1. Before the pipet controller is activated for the first time, the batteries should be charged.
2. If the pipet controller starts to indicate low battery level during work, connect it to the power supply to continue working.
3. Do not leave the pipet controller discharged for a long period of time.
4. In case of a planned longer break in using the device, it is recommended to remove the batteries from the battery compartment or charge them every few months.

9. Maintenance

Cleaning

The pipet controller does not require any maintenance. Its external parts may be cleaned with a swab moistened with isopropyl alcohol.

The nose piece and the pipet holder may be autoclaved at 121°C for 20 minutes. After autoclaving, dry the pipet holder. The filter included in the set may be sterilized by autoclaving at 121°C for not more than 15 minutes.

Ultra violet (UV) sterilization

The outer body of the pipet controller is UV resistant, which was confirmed by many tests. The recommended distance from the radiation source to exposed element should be not less than 50 cm.

Prolonged and very intense UV exposure can cause de-coloration of pipet controller parts, without affecting its performance.

Storage

The pipet controller should be stored in a dry place. The allowable storage temperature: -20°C to +50°C.

During breaks in the work the pipet controller can be stored on the wall hanger.

CAUTION:

Do not store the pipet controller with a filled pipet.

10. Components

The pipet controller set is supplied with the following components:

- Universal charger with set of adapters
- Membrane filter 0.2 µm (3)
- Batteries (3 AAA)
- Instruction manual
- Wall hanger
- Set of colored nose pieces (3: yellow, green, red)

11. Ordering information

The Corning® Stripettor™ Ultra Pipet Controller (Cat. No. 4099) comes with a universal charger and a set of adapters in different versions: EU, US, UK, and AUS. Choose your country's adapter and connect to the housing.

To mount the adapter, it should be inserted into the slots of the housing in the direction of the arrow (1), until you hear a click (Figure 6).

To remove or change the adapter, simply press the "PUSH" button in the direction of the arrow (2), hold the button down, remove the adapter in the direction of the arrow (3).

12. Spare parts for Corning Cat. No. 4099

Item in Figure 1	Description	Corning Cat. No.	Qty/Pk
F	Set of colored nose pieces (4: orange, yellow, green, red)	4997	4
G	Silicone pipet holder	4993	1
H	Membrane filter 0.2 µm	4991	5
	Membrane filter 0.45 µm	4992	5
L	Battery set - NiMH, AAA, 1.2V	4994	3
N	Universal power supply with set of adapters 9V: EU, US, UK, AUS	4996	1

13. Product disposal

According to Directive 2012/19/EU of the European Parliament and Council of 4th July 2012 on waste electrical and electronic equipment (WEEE) as amended, the Corning® Stripettor™ Ultra Pipet Controller is marked with the crossed-out wheeled bin and must not be disposed of with domestic waste.

Consequently, the buyer shall follow the instructions for reuse and recycling of waste electronic and electrical equipment (WEEE) provided with the products and available at the following link: www.corning.com/weee.

The crossed out wheeled bin symbol is printed in the product instruction manual and on the packaging to indicate that, in accordance with the requirements of Directive 2006/66/EC of 6th September 2006 on batteries and accumulators and waste batteries and accumulators, batteries must be disposed of accordance with national regulations.

For additional product information, visit www.corning.com/lifesciences, or call 1.800.492.1110. Customers outside the United States, call 1.978.442.2200 or contact your local support office.

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© 2016 Corning Incorporated

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476
[www.corning.com/
lifesciences](http://www.corning.com/lifesciences)

**Worldwide
Support Offices**

ASIA/PACIFIC

Australia/New Zealand

t 61 427286832

China

t 86 21 3338 4338

f 86 21 3338 4300

India

t 91 124 4604000

f 91 124 4604099

Japan

t 81 3-3586 1996

f 81 3-3586 1291

Korea

t 82 2-796-9500

f 82 2-796-9300

Singapore

t 65 6572-9740

f 65 6861-2913

Taiwan

t 886 2-2716-0338

f 886 2-2516-7500

EUROPE

France

t 0800 916 882

f 0800 918 636

Germany

t 0800 101 1153

f 0800 101 2427

The Netherlands

t 31 20 655 79 28

f 31 20 659 76 73

United Kingdom

t 0800 376 8660

f 0800 279 1117

**All Other European
Countries**

t 31 (0) 20 659 60 51

f 31 (0) 20 659 76 73

LATIN AMERICA

grupoLA@corning.com

Brasil

t (55-11) 3089-7400

f (55-11) 3167-0700

Mexico

t (52-81) 8158-8400

f (52-81) 8313-8589

For additional product information, visit www.corning.com/lifesciences, or call 1.800.492.1110.
Customers outside the United States, call 1.978.442.2200 or contact your local support office.

© Corning Incorporated. All rights reserved. Printed in Poland 4/16 CLS-EQ-113 REV2

Corning Incorporated, Daniszewska 4, 03-230 Warsaw
www.corning.com/lifesciences

Made in Poland